

Invia® Motion™ – Endure



Personal and convenient, offering clinical flexibility.



Invia® Motion™ — Endure NPWT Starter Kit

Advanced wound care facts:

- In the US, chronic wounds affect around 6.5 million patients a year and cost over \$25 billion to treat.¹
- The burden to treat wounds is growing rapidly due to increasing health care costs, an aging population and a sharp rise in the incidence of diabetes and obesity.¹
- Scheduling pick-up and returning equipment can cost an average of \$60 per home visit.²
- HCAHPS scores and patient readmission rates have direct impacts on reimbursement rates.³



Invia® Motion™ provides convenience

DME Providers

- Billable under HCPCS coding E2402
- Every patient receives their own starter kit containing all pump supplies
- Reduces need to retrieve, clean, and maintain equipment

Clinicians

- Every patient receives their own personal pump for the duration of therapy
- Visual and audible monitoring with step by step instructions
- Easy to use, only 4 buttons

Patients

- Size and weight promote mobility
- Quiet operation with quick canister replacement



¹Sen, Chandan K., Gayle M. Gordillo, Sashwati Roy, Robert Kirsner, Lynn Lambert, Thomas K. Hunt, Finn Gottrup, Geoffrey C. Gurtner, and Michael T. Longaker. "Human Skin Wounds: A Major and Snowballing Threat to Public Health and the Economy." *Wound Repair and Regeneration* 17.6 (2009): 763-71.

²Morrison Informatics, Inc. A Comprehensive Cost Analysis of Medicare Home Oxygen Therapy. A Study for the American Association for Homecare. June 27, 2006.

³<http://www.cms.gov/Medicare/Quality-Initiatives-Patient-Assessment-Instruments/HospitalQualityInits/HospitalHCAHPS.html> and <http://www.commonwealthfund.org/publications/blog/2013/feb/the-effect-of-medicare-readmissions-penalties-on-hospital>.

The Invia® Motion™ – Endure

Best suited for use in the home or other post-acute setting as it helps patients comfortably go about daily activities while promoting freedom and mobility during the entire length of treatment.



Features		
Pump category	Personal	
Duration of therapy	Customizable	
Pump weight	.88 lb	
Pump dimensions with canister	H x W x D 6.75 x 3.5 x 2.0 inches	
User interface	Soft membrane button touch	
Pressure options	-60, -80, -100, -125, -150, or -175mmHg	
Therapy flow	1 liter/min	
Therapy delivery mode	Constant or Intermittent	
Audible and visual patient notifications	Canister full, Leak in system, Battery low, Battery empty	
Canister capacity	150mL	
Power source	AC or battery	
Battery type	Internal NiMH	
Battery life	> 10 hours	
Battery recharge time	4 hours if empty	
Warranty	2 years	
ISO 9001 ISO 13485 CE (93/42/EEC), IIa	<p>Operation</p> <p>Transport/ Storage</p>	<p>Power supply Model: UE08WCP-090088SPA Input: 100-240VAC, 0.4A, 50-60 Hz Output: 9VDC, 0.88A Length of power supply cable: 67.71 inch</p>

Ordering information

Item #	Description	Case Quantity	HCPCS Code
0874004	Invia® Motion™ – Endure NPWT Starter Kit Includes – 1 Invia Motion – Endure pump – 1 150mL Canister/Tubing set – 1 Carrying case – 1 Power supply	5	E2402



Additional therapy supplies

Item #	Description	Case Quantity	HCPCS Code
0874001	150mL Canister/Tubing set with solidifier	10	A7000
0874004SB	Return Box	1	n/a
662151	Foam Dressing Kit - Small	5	A6550
662251	Foam Dressing Kit - Medium	5	A6550
662351	Foam Dressing Kit - Large	5	A6550
0877070	Invia® Gauze Kit – Standard ESI	15	A6550
0877061	Invia® Gauze Kit – Standard 10mm Flat Drain	15	A6550
0877063	Invia® Gauze Kit – Standard 15Fr Round Channel Drain	15	A6550
0877071	Invia® Gauze Kit – Large ESI	5	A6550
0877069	Invia® Gauze Kit – Large 10mm Flat Drain	5	A6550
0877074	Invia® Gauze Kit – Large 15Fr Round Channel Drain	5	A6550
0780025	External Suction Interface (ESI)	20	n/a

Replacement items

0874003	Carrying Case	10	n/a
0990243	Power Supply	1	n/a

Medical Vacuum Technology for Healthcare Professionals

Please contact us or your local Medela representative for details.

CAUTION: Federal US law restricts this device to sale by or on the order of a physician.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Medela Canada Inc.
4160 Sladeview Crescent, #8
Mississauga, Ontario
L5L 0A1
Canada
Phone +1 905 608 7272
Fax +1 905 608 8720
info@medela.ca
www.medela.ca

Medela Inc. USA
1101 Corporate Drive
McHenry, IL. 60050
USA
Phone +1 877 735 1626
Fax +1 815 307 8942
suction@medela.com
www.medelasuction.com

Local contact:

Medela and Invia are registered trademarks of Medela Holding AG.
Motion is a trademark of Medela Inc.