



EIKON™
illuminated
retractor system

Surgery. Redefined.



EIKON

illuminated retractor system

SUPERIOR INTRACAVITY VISUALIZATION

- BRIGHT, VOLUMETRIC AND UNIFORM ILLUMINATION
- THERMALLY COOL TECHNOLOGY

ADVANCED ERGONOMICS

- ACCOMMODATES MULTIPLE HAND POSITIONS
- HOLLOW CHANNEL PROVIDES CABLE MANAGEMENT

ATRAUMATIC BLADE

- SMOOTH EDGES REDUCE TISSUE TRAUMA

PRODUCT DESCRIPTION

Eikon Illuminated Retractor System

Set includes: B40XN, B40N, B50, B90, B90W, B135, B135W, B175, FC1S, ST1B

CATALOG NO.

B1S

Eikon Retractors

| Retractor | Dimensions | Catalog No. |
|-------------------|----------------|-------------|
| 40mm/Extra Narrow | (40mm x 12mm) | B40XN |
| 40mm/Narrow | (40mm x 20mm) | B40N |
| 50mm | (50mm x 30mm) | B50 |
| 90mm | (90mm x 30mm) | B90 |
| 90mm/Wide | (90mm x 60mm) | B90W |
| 135mm | (135mm x 30mm) | B135 |
| 135mm/Narrow | (135mm x 20mm) | B135N* |
| 135mm/Wide | (135mm x 60mm) | B135W |
| 175mm | (175mm x 30mm) | B175 |

Eikon Sterilization Tray


ST1B

WAVEGUIDES*

| | |
|---------------|--------|
| Wide/Flat | 104015 |
| Wide/Angled | 102521 |
| Narrow/Flat | 104008 |
| Narrow/Angled | 104001 |

ACCESSORIES

| | |
|---------------------------------------|------|
| Single Fiber Optic Cable, 10 ft | FC1S |
| Bifurcated Fiber Optic Cable, 11.5 ft | FC1B |



WITH INTELLIGENT PHOTONICS

www.invuity.com

Invuity, Inc. / 444 De Haro Street, San Francisco, CA 94107
Customer Service Tel: 866.711.7768 Fax: 888.614.8954

Invuity is a registered trademark, and Intelligent Photonics and Eikon are trademarks of Invuity, Inc. Patents www.invuity.com/patents. All rights reserved. ©2015 Invuity, Inc. LIT 10435 Rev F.

* Not included in B1S set

* Formerly known as Eigr Waveguide


intelligent photonics™