



EibergTM
Illuminated Retractor
System

Surgery. Redefined.

Eiberg™

Illuminated Retractor System

SUPERIOR INTRACAVITY VISUALIZATION

- BRIGHT, VOLUMETRIC AND UNIFORM ILLUMINATION
- THERMALLY COOL TECHNOLOGY

ADVANCED ERGONOMICS


- ACCOMMODATES MULTIPLE HAND POSITIONS COMFORTABLY
- HOLLOW HANDLE FACILITATES CABLE MANAGEMENT

ATRAUMATIC BLADE


- SUBTLE CURVATURES AND SMOOTH EDGES MINIMIZE TISSUE TRAUMA

UNIVERSAL TIP DESIGN

- HOHMANN-STYLE BLADE TIP FACILITATES SECURE PURCHASE INTO BONE



WITHOUT INTELLIGENT PHOTONICS™



WITH INTELLIGENT PHOTONICS™

PRODUCT DESCRIPTION

CATALOG NO.

Eiberg Retractors

Standard	(102mm x 27mm)	EIS
Standard Long	(140mm x 27mm)	EISL

WAVEGUIDES*

Wide/Flat	104015
Wide/Angled	102521

ACCESSORIES

Single Fiber Optic Cable, 10 ft	FC1S
Bifurcated Fiber Optic Cable, 11.5 ft	FC1B
Sterilization Tray, Large	110210

www.invuity.com

Invuity, Inc. / 444 De Haro Street, San Francisco, CA 94107
Customer Service Tel: 866.711.7768 Fax: 888.614.8954

Invuity is a registered trademark, and Intelligent Photonics and Eiberg are trademarks of Invuity, Inc.
Patents www.invuity.com/patents. All rights reserved. ©2015 Invuity, Inc. LIT 12090 Rev B

INVUITY®
intelligent photonics™

CE 0086